Acryl BH

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Instructions for use

Application:

Fabrication of partial and full dentures using the flasking and pressing technique

Characteristics:

- Colour stable
- Easy to use

Mixing ratio:

2.5 parts powder by weight

1 part liquid by weight

Approx. **25 g powder** and **10 g liquid** are required for a full denture.

Technique:

- Coat the plaster surfaces with alginate separating agent.
- After adding the powder to the liquid, spatulate Acryl BH thoroughly to a homogeneous mixture in a mixing vessel for 30 secs.
- Allow a doughing time of approx. 15 minutes. The acrylic dough then has a working time of approx. 30 minutes.
- At the end of the doughing time knead the acrylic dough thoroughly for approx. 30 seconds.
- When the Acryl BH becomes viscoelastic, pack it into the prepared flask and press.
- Hold the pressure for approx. 15 minutes.
- Secure the flask in a clamp for polymerisation.
- The denture teeth should be roughened, provided with undercuts and coated with monomer to ensure they are securely retained.

Polymerisation:

Polymerize Acryl BH in the closed flask with the clamp for 20 minutes in boiling water. Allow the flask to cool to room temperature in the clamp and then deflask in the usual way.

Repairs:

Acryl BC or Acryl S can be used for repairs. Before applying the acrylic dough, the fractured surfaces should be roughened, cleaned and coated with monomer. Polymerize in a pressure vessel for 15 minutes at a pressure of 2 - 5 bar with a water temperature of 45 °C.

Storage:

Store at 10 - 25 °C. Do not use after the expiry date.

Important:

Should the room temperature deviate from 23 °C: higher temperatures reduce the working time and lower temperatures extend the working time.

Caution:

Avoid skin contact with unpolymerized material and inhalation of monomer fumes.

Opened containers: Do not use up after the expiry date.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare.

Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-Indications/interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Disposal:

Powder can be disposed of as domestic waste. Liquid should be disposed of as hazardous waste.

Troubleshooting

Problem	Cause	Solution
acrylic does not cure completely	<ul style="list-style-type: none">– incorrect mixing ratio	<ul style="list-style-type: none">– adhere to instructions– adhere to mixing ratio
acrylic shows a whitish tint at the bottom	<ul style="list-style-type: none">– incorrect or insufficient insulation– model was not soaked in water	<ul style="list-style-type: none">– insulate properly– check insulation agent, if necessary– soak model in lukewarm water for ca. 10 min.
bubbles form at the bottom of the acrylic	<ul style="list-style-type: none">– model was not soaked in water	<ul style="list-style-type: none">– soak model in lukewarm water for ca. 10 min.
surface of the acrylic appears white and porous	<ul style="list-style-type: none">– acrylic was injected too late respectively too much time passed until it was set into the pressure pot	<ul style="list-style-type: none">– adhere to the processing times
color deviations between repairs and original restoration	<ul style="list-style-type: none">– incorrect mixing ratio	<ul style="list-style-type: none">– adhere to mixing ratio
acrylic breaks when lifted	<ul style="list-style-type: none">– insufficiently insulated	<ul style="list-style-type: none">– check insulation
acrylic / restoration shows whitish streaks	<ul style="list-style-type: none">– badly mixed / too much powder used	<ul style="list-style-type: none">– carefully mix acrylic dough with a spatula, adhere to mixing ratio

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