# Acryl N

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## Directions for Use

### Indications:

Adding the saddles to chrome cobalt denture bases. Repairing, extending and relining dentures using liquid resin.

### Properties:

- Colour stability
- High mechanical values

### Advantages:

Easily processed.

### Mixing ratio:

- 10 parts by weight powder
- 7 parts by weight liquid

### Processing:

- Apply separating agent to the plaster.
- Add the powder to the liquid, allow 15 seconds for it to soak in completely and spatulate the Acryl N thoroughly in a mixing bowl.
- Initiates fluid state, Acryl N can be poured into, for example, an overcast. Once the plastic phase sets in, the Acryl N will no longer flow out of the overcast and can be shaped.
- To prevent acrylic teeth breaking off the denture the teeth, as is standard procedure for cold-curing resins, must be roughened, provided with undercuts and wetted with monomer or Cross-Liquid.
- After a maximum of 6 minutes (at a room temperature of 23°C), must be placed in a pressurized polymerization unit for polymerizing.

### Polymerization:

Acryl N should be polymerized in a pressurized polymerization unit for 15 minutes, at a pressure of 2-5 bars and a water temperature of 45°C.

### Repairs - Please note:

Repairs can be carried out with Acryl S, Acryl N and Acryl LT.

Prior to applying the resin dough, the fractured surfaces must be roughened, cleaned and wetted with monomer. The resin should be polymerized in a pressurized polymerization unit for 15 minutes, at a pressure of 2-5 bars and a water temperature of 45°C.

### Storage:

Keep at 10 - 25°C.

### Important:

Should the room temperature deviate from 23°C: Higher temperatures reduce and lower temperatures prolong the working times.

### Hazard note:

Avoid prolonged skin contact with unpolymerized material and inhalation of monomer vapours.

Once the container has been opened: Do not use after the expiry date.

### Disposal:

The powder can be disposed of with household waste. Dispose of the liquid in a special refuse tip.

### Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

### Contra-indications / interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

### Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

## Troubleshooting

Problem	Cause	Remedy
Acrylic does not fully cure	Incorrect mixing ratio	<ul style="list-style-type: none"> <li>- Adhere to the instructions for use</li> <li>- Adhere to the mixing ratio</li> </ul>
Acrylic is whitish on the fitting surface	<ul style="list-style-type: none"> <li>- incorrect or poor application of separating agent</li> <li>- model not soaked</li> </ul>	<ul style="list-style-type: none"> <li>- Apply ample separating agent</li> <li>- Check the separating agent if necessary</li> <li>- Soak the model for approx. 10 min. in hand-hot water</li> </ul>
Acrylic form bubbles on the fitting surface	Model not soaked	<ul style="list-style-type: none"> <li>- Soak the model for approx. 10 min. in hand-hot water</li> </ul>
Acrylic surface white / porous	Recommended time for injecting acrylic exceeded or delay in inserting into the pressure vessel	Adhere to the processing times
Colour differences with repairs	Incorrect mixing ratio	Adhere to the mixing ratio
Acrylic breaks when lifted	Poorly separated	Check the separating agent
Acrylic / denture has whitish streaks	Incorrectly mixed / too much powder used	Thoroughly spatulate the acrylic dough / adhere to the mixing ratio
Acrylic sets too fast	Temperatures too high Working stages too long	Observe temperatures / working times Storage: 10 - 25°C

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