# DIE

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## Directions for use

### Range of Applications

DIE is excellent for fabricating precision models for all types of crown and bridge-work as well as for dimensionally stable reproduction of shoulders and partial crown preparations.

**DIE Separator must be applied when casting polyether impression, resp., duplicating materials.** When using hydrocolloid impression materials or duplicating materials, the moulds should be dried with compressed air.

**Do not use with alginates.**

### Properties

DIE is flowable, accurate in detail and dimensionally stable. It has high edge strength and abrasion resistance and is easy to trim and saw.

### Comprises

DIE is supplied in glass bottles (200 ml base material, 100 ml hardener), together with 400 g filler.

### Mixing DIE

Briefly shake the bottles containing both the base material and hardener prior to mixing. Note the coloured marking on the cap. Close the glass bottles carefully immediately after use. **Do not allow DIE to come into contact with water** (foam).

1. Fill the required amount of base material in the dispensing and mixing container supplied. Then add the correct quantity of

filler and thoroughly spatulate the mixture. Add the correct quantity of hardener and spatulate the mixture again thoroughly.

### Mixing ratio:

Base: hardener = 2 : 1.

Approx. 15 g DIE and filler is required for a full dental arch.

10 g base

1 level measuring scoop of filler (approx. 14 g)

5 g hardener

2. After mixing for approx. 30 sec., pour the DIE in a thin stream into the cleaned and dried impression. The material remains flowable for approx. 2 min. at 20°C.

3. In order to prevent incorporation of air voids, DIE may be hardened for 15 mins., after pouring, in a dry pressure vessel at 2-4 bars.

DIE is sufficiently hard after 1 hour that the arch can be trimmed, drilled and ground.

### CAUTION!

Do not overheat DIE. Steam-clean DIE (max. 15secs.).

### Disposal

Dispose of in a special refuse tip, in accordance with local regulatory requirements.

Carefully carried out working procedures ensure that the die meets the most stringent requirements.

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