Nano Composite

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Instructions

Nano Composite is a light-curing nano-composite for the adhesive filling technique. It contains an ultrafine, radiopaque glass filler. The composition with ultrafine filler leads to exceptionally homogenous restorations. An optimal adaption of the color is the result of an optimally adjusted chameleon effect. The guidelines and requirements of DIN EN ISO 4049 apply. Nano Composite should be used in connection with the system components Etching Gel, Bond LC and Bond LC Activator.

Nano Composite is available in syringes and compules. The compules are for single use only. Please do not reuse them, as this makes it impossible to rule out contamination and germ formation.

Composition

Monomer matrix

Durethane dimethacrylate, butanedioldimethacrylate, isopropylidene-bis [2(3)-hydroxy-3(2)-(4-phenoxy)-propyl]bismethacrylate [Bis-GMA]

Filler content 83.5 % by weight (66.5 % by volume)
Anorganic filler (28 nm -15 µm)

Indications

- Direct anterior and posterior tooth restorations in Black’s classes I, II, III, IV and V.
- Indirect restorations such as inlays, onlays and veneers
- Extended fissure sealing on molars and premolars
- Building up stumps
- Splinting of loosened teeth
- Corrections of shape and color to enhance aesthetics

Available colors

Incisal masses:
InW (white), InTr (transparent), InBI (bleach), InU (universal), IR (red)

Dentin masses: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, D3

Opaque dentin: OD A2, OD B2

Selection of colors

The selection of the colors depends on the position and the size of the cavity as well as using the color of the remaining tooth substance and/or the color of the remaining dentition.

The color selection for the cervical sectors should be made among the different darker yellow and red shades.

For the dentin sector, select colors from among the yellow, grey or red shades.

For the incisal sector, use material in lighter or transparent shades.

Application

Preparation of the natural tooth
Before starting the treatment, clean the hard tissue with fluoride-free polishing paste. In order to avoid a possible pulp reaction, always prepare a liner in cases with exposed dentin (e.g. compound containing calcium hydroxide).

Cavity preparation
Gently prepare the hard tissue according to the rules of the adhesive technique. When working on anterior teeth, bevel all enamel edges. Do not bevel the edges when working on posterior teeth. Avoid feathered margins. Next, clean the cavity with water spray, remove all residue and dry. Complete drying is necessary. We recommend the use of a kofferdam.

Pulp protection / liner
When using an enamel-dentin adhesive the liner can be foregone. In case of very deep cavities near the pulp, line with a calcium hydroxide compound accordingy.

Design of approximal contacts
In cavities with approximal sections, set in a transparent matrix and fix.

Etching and bonding

Apply Etching Gel to the preparation surface, etch enamel for 30 sec. and/or etch dentin for 15 sec. Thoroughly rinse off etching gel with water. Avoid drying out the dentin. Dried, etched enamel surfaces have chalky-white appearance and must not be contaminated before the bonding is applied. If the preparation area is contaminated by saliva, rinse and dry again. Etch once again, if necessary.

Shake Bond LC well before use. Liberally apply Bond LC to the dentin and enamel surfaces with a brush. Brush in intensively for about 30 sec. Afterwards, carefully dry with oil-free compressed air for about 15 sec. Cure the Bond LC layer for 20 sec. with a dental polymerization lamp, before applying a second layer of Bond LC in the same way.

Applikation of the composite

Take the required amount of composite out of the application tip. Fill the cavity with the material and shape as needed, using conventional metal instruments. A layer thickness of 2 mm must not be exceeded. Due to the oxygen in the ambient air, a thin dispersion layer will remain on top of each layer. This dispersion layer forms the chemical bond between the layers and must not be touched or contaminated with moisture.

Polymerization
The polymerization time is 40 sec. per layer for colors with a commercial polymerization system, an LED polymerization lamp or 2 x 3 sec. with a plasma polymerization system. Hold the light-guide as close to the surface of the filling as possible. Fillings with several surfaces should be polymerized from the direction of each surface.

Finishing
Nano Composite can be finished and polished immediately after polymerization. Suitable are finishing diamonds, flexible separating discs, silicone polishers and polishing brushes. Check occlusion and functional movements, correct if necessary. Finally, polish with suitable polishing pastes.

Indirect method
Cavity preparation:
A tooth-conserving preparation with only little divergence in the cavity walls should be preferred. All inner edges and angles must be rounded, avoid feather margins. Design a planar cervical shoulder, do not bevel it. Block out any unavoidable undercut with glass ionomer cement. For the preparation, use slightly cone-shaped diamond grinders with rounded edges. Cover dentin close to the pulp with a thin layer of a calcium hydroxide compound. Any liner materials that contain eugenol are contra-indicated.

Impression and temporary restoration
After taking the impression, a temporary resin restoration is produced. Fix temporaries only with eugenol-free cement.

Production of inlays, onlays and veneers

At the lab, pour an extra-hard plaster into the impression. Once the die has hardened, remove the impression from the die. Block out undercuts and insulate the die with an oil-free insulator. Build up the inlay on the die layer by layer. First, build up approximal and deep occlusal areas. Each layer should have a maximum thickness of 2 mm.

Polymerize with a commercial polymerization system (e.g. HiLite Power, Heraeus Kulzer, intermediate polymerization 90 sec., final polymerization 180 sec.). Lift off the completed restoration from the die, finish and polish to a high gloss. Clean the restoration thoroughly with water and soap. Rinse with air/ water spray and dry.

Inserting inlays, onlays and veneers

Remove the temporary restoration and clean the cavity. Set in a kofferdam, then clean and dry the surface of the tooth. Check the restoration’s fit with slight pressure. Avoid rough handling. If necessary, improve the fit by grinding the inside surface. Do not check the occlusion in this first fitting, as the restoration could break. Apply Etching Gel onto the preparation surface. Etch enamel for 30 sec. and/or dentin for 15 sec. Thoroughly rinse off the etching gel with water. Avoid drying out the dentin. Dried, etched enamel surfaces have chalky-white appearance and must not be contaminated before the bonding application. If the preparation area is contaminated with saliva, rinse and dry again. Etch once again, if necessary.

Apply a thin layer of Bond LC bonding agent to the etched areas and rub in. Polymerize with a dental polymerization lamp for 40 sec. During polymerization, a dispersion layer is formed which must not be removed, as it forms the chemical bond to the filling material.

Fixing the restoration

The restoration is fixed with a commercially available, dual-curing fixing composite. Please adhere to the manufacturer’s instructions.

Special instructions

- The working range under a surgical light is 2 min.
- In case of extensive restorations, the surgical light should be moved away from the working area temporarily to avoid premature curing of the composite. Alternatively, the material can be covered with a light-tight foil.
- Use a light polymerization system with an emission range of 350 -500 nm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer’s instructions.

Storage

Store at temperatures between 10 °C and 25 °C (50 °F to 77 °F) . Avoid direct sunlight. Screw the cap back onto the syringe tightly after each use. Let the material reach room temperature before use. Withdraw the plunger slightly after use to keep the outlet from becoming plugged. Do not use after the expiration date (see label on the syringe). Only for use in dentistry. Keep out of children’s reach. This product has been developed for the specific use

illustrated above. Only process as described in these instructions. The manufacturer will not be held liable for any damages that result from improper use or improper processing.

We recommend to forego the use of syringe heaters.

Side-effects:
With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e.g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. In order to avoid a possible pulp reaction, always prepare a liner in cases with exposed dentin (e.g. compound containing calcium hydroxide).

Contraindications / interactions:

If a patient has known allergies against or hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient’s mouth before using the product. Phenolic substances (e.g. eugenol) inhibit polymerization. Therefore, these materials (e.g. zinc oxide eugenol cements) must not be used as liners.

* Vita is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.

Hazard statement

May cause an allergic skin reaction.

Precautionary statement

Use protective gloves/ protective clothing /eye protection /face protection. If skin irritation or rash occurs: Get medical advice/attention.

Trouble shooting		
Trouble	Cause	Remedy
Composite does not cure.	Light output of the polymerization lamp is insufficient.	Check the light output and, if necessary, replace the light source.
	Emitted wavelength range of the polymerization lamp is insufficient.	Consult the manufacturer of the polymerization lamp. Recommend wavelength range 350 - 500 nm.
Composite inside the syringe is sticky and soft; clear liquid separates from the material inside the syringe.	Material was stored at temperatures ≥ 25 °C (≥ 77 °F) for a longer period of time.	Please note the storage temperatures: Store at temperatures between 10 °C (50 °F) and 25 °C (77 °F). Short-duration storage in a refrigerator.
	Material was stored in the syringe heater for too long.	Never leave a syringe longer than one hour per use in a syringe heater.
Composite appears to be compact and hard inside the syringe.	Material was stored at temperatures < 10 °C (< 50 °F) for a longer period of time.	Let composit reach room temperature before use. If necessary, use a syringe heater.
	Syringe was not closed correctly, composite has started to polymerize.	Each time after removing material from the syringe, correctly cap and close the syringe.
Inlay/onlay does not stay in place after insertion.	Restoration appears too opaque to be fixed with a purely light-polymerizing composite.	Use dual-polymerizing fixing composite.
Composite does not cure all the way through (dark or opaque colors).	Layers per polymerization cycle were too thick.	Do not exceed max. layer thickness of 2.0 mm per layer.
Restoration has a yellowish tint when compared to the color reference.	Insufficient polymerization of the composite layers.	Repeat polymerization cycle several times; at least 40 sec.

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