Root Cem Viscous

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Instructions

Description
Root Cem Viscous is a dual-curing fixing composite. The material is fluorescent and radiopaque.

Indications

Root Cem Viscous is indicated for:

- the adhesive fixation of indirect restorations such as veneers, inlays, onlays, crowns and bridges. The material is suited to fixing restorations made from glass ceramics, zirconium dioxide, composite, zirconia-reinforced composite, high-precious metal, non-precious metal, titanium.
- the adhesive fixation of root pins.
- building up stumps (as core build up material).

Contraindication

The use of **Root Cem Viscous** is contraindicated

- if a patient has known allergies against one or several of the product’s components
- if it is impossible to clean and dry the treatment site as necessary or if the prescribed technique cannot be used.

Processing range

Processing and curing times are dependent on the ambient temperature. As soon as **Root Cem Viscous** has been removed from the automix cartridge, the following times are to be considered:

	Ambient temperature ca. 21 °C (ca. 69 °F)	Intraoral 37 °C (98,6 °F)
Processing range	3- 4 min.	ca. 2 min.
Curing time incl. processing time	ca. 7- 9 min.	ca. 3- 4 min.

Light polymerization

Light intensity	Optical wavelength range	Exposure time per surface
> 500 mW/cm²	350- 500 nm	20 sec.

Note: After excessive material has been removed, please polymerize all edges again for 20 sec.

Mixing ratio

With the automix cartridge and the mixing tip, **Root Cem Viscous** is always available in the optimal mixing ratio of 1:1.

Composition

Glass powder, diurethane dimethacrylate, prepolymer, tetramethylene dimethacrylate, silicon dioxide, dibenzoyl peroxide

Filter content:

75 % by weight of inorganic filler (0.005 - 5.0 µm)

Application advice

- The material should have reached room temperature before processing, as the material can be harder to press out and mix at refrigerator temperatures.
- Process **Root Cem Viscous** immediately after pressing it out of the automix cartridge and set in the restoration quickly.

Use

Removal of the temporary restoration
Remove the temporary restoration, clean cavity or stump, remove any residue of the cement that has fixed the temporary restoration and rinse well with water. Subsequently, dry with oil-free compressed air.

Try-in of the restoration

- Try in the restoration, check fit and contact points.
- Clean the surfaces to be glued thoroughly with alcohol (e. g. isopropanol). Subsequently, dry with oil-free compressed air.
- Please note:** PMMA materials must **NOT** be cleaned with alcohol. Clean such objects e. g. with a disinfectant based upon quaternary ammonium salts.
- Moisture ingress and contamination must be avoided.
- In case of delicate ceramic restorations, the occlusion test should be carried out very carefully, as before the final fixation of the restoration, there is always a danger of fracturing it.
- If necessary, corrections are carried out with fine diamonds at a medium revolution speed. Any surfaces that were ground must be polished again.
- It is essential for the adhesive fixation of a restoration to dry the working area well (preferably with cotton rolls or parotid gland patches). We recommend a complete isolation with e.g. a kofferdam.

A) Fixation of indirect restorations

A1) Pretreatment of the natural tooth

- Apply the etching gel to the preparation surface. Etch enamel for 30 sec. and/or dentin for 15 sec.
- Afterwards, rinse off etching gel well with water.
- Wet the surface with a mixtu of Bond LC and Bond LC Activator (1:1). Alternatively, a different adhesive system can be used. Please adhere to the manufacturer’s instructions.
- We recommend to apply the adhesive with a microbrush.
- Blow off water and solvent residue with a gentle airstream.
- The mixture of Bond LC and Bond LC Activator (1:1) is applied again in a second layer and blown off.

A2) Pretreatment of the restoration

- Generally follow the instructions of the restoration material’s manufacturer.

Pretreatment of	Procedure
<ul style="list-style-type: none">Metals	<ul style="list-style-type: none">Sandblast the inside surfaces of the restoration with 30- 50 µm aluminium oxide at a maximum pressure of 2 bar Clean with alcohol and dry the restoration Apply Universal bonder evenly with a thin brush, afterwards let it dry for ca. 1 min. Note: For an optimal bond, do not clean the surfaces with phosphoric acid!!!
<ul style="list-style-type: none">Lithium disilicate glass ceramics	<ul style="list-style-type: none">Etch with 5% hydrofluoric acid for 20 sec. or in accordance to the glass ceramics manufacturer’s instructions. Thoroughly rinse the restoration with water and dry with oil-free air. Apply Universal bonder evenly with a thin brush, afterwards let it dry for ca. 1 min.
<ul style="list-style-type: none">Composite	<ul style="list-style-type: none">Sandblast the inside surfaces of the restoration with 30- 50 µm aluminium oxide at a maximum pressure of 2 bar Clean with alcohol and dry the restoration Apply Bond LC evenly with a thin brush, blow-dry with a gentle airstream, but do not polymerize.
<ul style="list-style-type: none">PMMA materials	<ul style="list-style-type: none">Sandblast the inside surfaces of the restoration with 30- 50 µm aluminium oxide at a maximum pressure of 2 bar Clean with a suitable agent, e. g. disinfectant based on quaternary ammonium salts and dry the restorations Apply Universal bonder evenly with a thin brush, afterwards let it dry for ca. 1 min. Note: Avoid contact with alcohol as it can lead to cracks.

A3) Application of Root Cem Viscous

- Press out a small amount of paste and discard. This is done to bring in line the base paste and catalyst inside the dual cartridge.
- Set a new disposable automixing tip onto the dual cartridge. Please do not use the first bit of mixture coming out of the mixing tip for fixation.
- Press out **Root Cem Viscous** from the automix cartridge and apply the desired amount directly onto the restoration. Spread evenly.
- As the fixing material will cure inside the mixing tip, it can be used as a cap for the automix cartridge until its next use to protect the contents.

A4) Inserting the restoration and removal of excess material

- Bring restoration into position with slight pressure. Fix and hold.
- Remove any excess fixing composite immediately with a disposable brush, a cotton pellet, dental floss or an implant scaler. Pay special attention to the immediate removal of excess material in hard-to-reach areas (approximal, gingiva margin).
- In case of composite or glass ceramic restoration, light-polymerize for 20 sec. from each tooth surface after removal of excess material. Chemical curing is finished after 6 min.
- Please also light-polymerize all edges for 20 sec. after removal of excess material.
- To avoid the formation of an oxygen inhibition layer, the margins of the restoration can be covered with a layer of glycerin gel/air

blocker immediately afterwards. This layer can be rinsed off with water after the polymerization has been completed.

A5) Finishing the restoration

- If applicable, remove the retraction thread.
- Check occlusion and functional movements, correct if necessary.
- Polish margins and restoration surfaces.

B) Fixation of root pins

- Prepare the root canal and cavity as usual.
- Clean the root canal with a 3% hydrogen peroxide solution (H₂O₂) or with a 2.5- 5.25% sodium hypochlorite solution (NaOCl).
- Immediately rinse with water and dry with paper tips.
- Select an appropriate root pin and try in.
- Clean the root pin with alcohol.
- Moisten the root pin with mixed Bond LC and Bond LC Activator (1:1).
- Etch cavity (coronal and radicular part of the preparation) with Bond LC (enamel 30 sec., root dentin 15-30 sec.).
- Thoroughly rinse off etching gel with water, remove all acid residue and dry the cavity with paper tips. Do not blow dry with air in order to keep the dentin moist.
- Rub a mixture of 1:1 Bond LC and Bond LC Activator with a thin microbrush onto the occlusal surface of the stump inside the root canal for 15 sec. Remove any excess material with paper tips.
- Attention:** Ensure that the micro brush reaches all areas of the root canal and that the bonding agent is spread evenly. The grip of the brush should not touch the tooth surface or get jammed. Blow off water and solvent residues with air.
- Position root pin to check the depth of the canal and to press the adhesive further into the dentin tubuli. Afterwards, remove the pin.
- Set a new disposable automixing tip onto the dual cartridge. Please do not use the first bit of mixture coming out of the mixing tip for fixation. Fill the canal from the bottom of the cavity to the corona with Root Cem Viscous. Do not use lentils or other rotating instruments.
- Coat the surface of the pin with a sufficient quantity of Root Cem Viscous and place the pin completely into the canal. Let any excess fixing composite drain. Fix the pin and remove surplus cement carefully.
- After setting it in, light polymerize for 40 sec. Afterwards, start the restoration.

C) Building up a stump

- If the coronal surface has not yet been etched as described under B, apply Etching Gel to the dentin and etch for 15 sec.
- Thoroughly rinse off the etching gel with water, remove all acid residue.
- Extract water and carefully dry the surfaces with paper tips.
- Do not blow dry with air to keep the dentin moist.
- We recommend to apply the adhesive with a microbrush.
- Blow off water and solvent residue with a gentle airstream.
- Apply Bond LC without activator and light polymerize for 20 sec.
- Set a new disposable automixing tip onto the dual cartridge and apply an appropriate amount of **Root Cem Viscous** onto the root pin.
- If a stump matrix is used, fill it and set it over the pin.
- Remove excess material and light polymerize for 40 sec.
- Chemical curing is complete after 6 min.
- Finish and polish. An impression can be taken immediately afterwards.

Warning

Contains tetramethylene dimethacrylate, dibenzoyl peroxide
Attention: May cause an allergic skin reaction. Avoid breathing vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

Contraindications:

If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient’s mouth before using the product. Unpolymerized resin may lead to allergic skin reactions. Therefore, the user should consider protective measures. In case that irritations occur or in case of a known allergy against one of the components listed above, please refrain from using the product.

Interactions with other substances

Phenolic substances (like e.g. eugenol) inhibit the polymerization. Therefore, do not use any liners that contain such substances (e. g. zinc oxide-eugenol-cements). Oxidative disinfectants (e. g. hydrogen peroxide) can interact with the polymerization system and thus negatively influence the polymerization. Therefore, do not use an oxidative disinfectant on the automix cartridge. The automix cartridge can be disinfected by e.g. wiping down with medical alcohol.

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