

Bulk-Fill Composite

Light-curing, radiopaque posterior composite for easy and fast handling in bulk-fill technique.

The material is both flowable and highly viscous.

Features	Benefits
Up to 4 mm layers	Highly efficient handling
Universal dentin shade with natural colour adjustment	Saves time
Optimised flow behaviour	Excellent adaption to cavity bottom and cavity walls
Excellent physical properties	High quality, long-lasting posterior restorations
High Radiopacity	Ideal for x-ray checks
Doesn't contain Bis-GMA	Optimised biocompatibility



Shade: Universal Dentin
Packaging: 2 g syringe or 0.25 g Filli
Accessories: Application tips

Technical Data:

Vickers hardness 737 MPa
Flexural strength 134 Mpa
E-Module 11090 MPa
Compressive strength 381 MPa
Filler content Approx.77%
Radiopacity 210% Al

Version 10/2014