


Gingiva Composite PASTE & FLOW

Light-curing Gingiva Composite in two consistencies (paste and flow) for gingival restorations and reproductions, implant suprastructures, Crowns & Bridges as well as characterization of full & partial dentures. The material comes in four different shades for natural looking gingival reproductions.

Features	Benefits
4 shades in paste and flow	Natural looking reconstructions of gingival parts
Abrasion-resistant	For highly stable restorations
Excellent polishability	Low plaque affinity
Compatible with light-curing veneering composites (see C&B Composite Paste & Flow)	Perfect addition to dentine shades



Shades:
Light – Medium – Dark – Super Dark

Packaging:
5 g syringe PASTE
3 g syringe FLOW
(accessories for flowable:
10 application tips recommended)

Storage: 10-25°C
Shelf-life: 4 years



Technical data:	Paste	Flow
Hardness	400 Mpa	370 MPa
Flexural strength	95 MPa	140 MPa
E-Module	6900 MPa	6500 MPa

Version 10/2014