

## Polish 1 + Polish 3

Premium quality and fine grained diamond polishing pastes for professional final polishing of aesthetic direct and indirect composite and ceramic restorations, veneers, inlays etc.

Features	Benefits
High diamond content	Quick polishing results with a minimum amount of paste
Fine grained particle size 1 und 3 $\mu\text{m}$ , respectively	High-gloss shine
Firm consistency	Easy to use
Universal use	For dentist an dental lab



Packaging: 2 g syringe or 10 g jar

Storage: 10-25°C

NB: To be applied with a goat hair brush

Version 05/2016