

Bond LC

Light-curing, one-component bonding suitable for enamel as well as dentine. Suitable as a primer for ceramic and composite inlays as onlays, veneers, crowns and bridges well as for root canals before filling or fixing root pins.

Features	Benefits
Primer and Bond in one bottle	Quicker and more reliable to use
Ethanol based	Very stable
Optimal viscosity	Good biocompatibility
Hydrophilic	Easy to apply. Penetrates the dentine tubule well



Packaging size: 5 ml bottle
Storage: 10-25°C

Bond LC Activator

To be mixed with Bond LC: Dual-curing, two-component bonding suitable for dentine as well as enamel. Suitable as a primer for inlays as well as for root canals before fixing or filling root pins.

Features	Benefits
Primer and Bond in one bottle	Quicker and more reliable to use
Ethanol based	Good biocompatibility
Optimal viscosity	Very stable
Hydrophilic	Readily wettable
Dual-curing	Safe curing without light



Packaging size: 5 ml bottle
Storage: 10-25°C
Accessories: Dental brush holder, mixing pallet, dental brushes

Technical Data:

Bonding stability to dentin and etched enamel: 30 MPa
Bonding stability to non-precious metals (Co-Cr): 24 MPa

Version 04/2016